

CONNECTIVITY OPTIONS FOR REXTON HEARING AIDS

	Audio streaming				Hands-free calls via Headset Mode			Control via Rexton App
	Bluetooth Classic	Bluetooth LE Audio	ASHA	MFi	Bluetooth Classic	Bluetooth LE Audio	MFi	
Reach R Plus	■	-	-	-	■	-	-	■
Reach R -Li / R-Li T	-	■	■	■	-	■	■	■
Reach SR	-	■	■	■	-	■	■	■
Reach B-Li P / B-Li HP	-	■	■	■	-	■	■	■
Reach B-Li Rugged	-	■	■	■	-	■	■	■
Reach InoX-CIC Li	-	-	-	-	-	-	-	■
Reach Custom CIC-Li	-	-	-	-	-	-	-	■
Reach Custom	-	-	-	-	-	-	-	■
BiCore R-Li / R-Li T	-	-	■	■	-	-	■	■
BiCore R312	-	-	■	■	-	-	■	■
BiCore SR	-	-	■	■	-	-	■	■
BiCore B-Li Rugged	-	-	■	■	-	-	■	■
BiCore B	-	-	■	■	-	-	-	■
BiCore Custom Li	-	-	■	■	-	-	■	■
ReCharge R-Li	-	-	-	-	-	-	-	■
ReCharge B-Li	-	-	-	-	-	-	-	■