

REXTON HEARING AID COMPATIBILITY WITH CONNECTIVITY STANDARDS

Apple devices	Android devices	all Bluetooth devices
Bluetooth Classic	Bluetooth Classic	Bluetooth Classic
MFI	ASHA	
	LE Audio	

	Audio Streaming				Hands-free calls via Headset Mode			Control via Rexton App
	LE Audio	MFI	ASHA	Bluetooth Classic	LE Audio	MFI	Bluetooth Classic	
Reach R Plus				■			■	■
Reach R -Li / R-Li T	■	■	■		■	■		■
Reach SR	■	■	■		■	■		■
Reach B-Li Rugged	■	■	■		■	■		■
Reach InoX-CIC Li								■
Reach Custom								■
BiCore R-Li / R-Li T		■	■			■		■
BiCore R312		■	■			■		■
BiCore SR		■	■			■		■
BiCore B-Li Rugged		■	■			■		■
BiCore B		■	■					■
BiCore Custom Li		■	■			■		■
ReCharge R-Li								■
ReCharge B-Li								■