

# REACH CUSTOM CIC-LI


## CUSTOMIZED INVISIBLE RECHARGEABILITY



- > Our latest rechargeable custom device for a perfect fit in the ear canal
- > With Voice Stabilizer the contrast between speech and noise in every conversation is controlled for easier speech intelligibility
- > OneMic Directionality enables the benefits of a directional microphone
- > A discreet design with easy handling makes it virtually invisible
- > Portable charger with built-in power bank
- > Rechargeable solution for a full days use:
  - > Custom CIC-Li delivers up to 35h runtime
- > Lifeproof to withstand everyday challenges

- Faceplate colors
- Black
  - Mocha
  - Dark Brown
- Shell colors
- Black
  - Mocha
  - Dark Brown

Fitting range is based on 2 ccm coupler:



- |  |  |  |
|--|--|--|
| <ul style="list-style-type: none"> <li>• 50 dB / 113 dB SPL (2 ccm coupler)</li> <li>• 60 dB / 124 dB SPL (Ear simulator)</li> </ul> | <ul style="list-style-type: none"> <li>• 55 dB / 118 dB SPL (2 ccm coupler)</li> <li>• 65 dB / 127 dB SPL (Ear simulator)</li> </ul> | <ul style="list-style-type: none"> <li>• 60 dB / 124 dB SPL (2 ccm coupler)</li> <li>• 70 dB / 133 dB SPL (Ear simulator)</li> </ul> |
|--|--|--|

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	36/18/6	32/16/6	24/12/4	16/8/4
LED	-	-	-	-	-
T-coil	-	-	-	-	-
Bluetooth	-	-	-	-	-
Tinnitus Sound function	•	•	•	•	-
Tinnitus Notch function	•	•	•	•	-
Wireless sync	•	•	•	•	•
IP68	•	•	•	•	•
XPhone	•	•	•	•	•
<b>Technology Highlights</b>					
Multi-Voice Focus	-	-	-	-	-
Focus Beams	-	-	-	-	-
Voice Stabilizer	•	•	•	•	•
Speech Preservation Technology	-	-	-	-	-
My Voice 2.0	-	-	-	-	-
Motion sensor for scene analysis	-	-	-	-	-
<b>Directionality</b>					
Auto Dir	•	•	•	-	-
Directional Adaptive	•	•	•	-	-
iOmni Sound Locator	-	-	-	-	-
Auto iFocus 360	-	-	-	-	-
iFocus 360	-	-	-	-	-
Auto Focus 360	-	-	-	-	-
Auto Stereo iLock	-	-	-	-	-
OneMic Directionality	•	•	•	-	-
<b>Noise Reduction</b>					
Noise Management	•	•	•	•	•
Sound Smoothing	•	•	•	•	•
Wind Noise Reduction	•	•	•	•	•
Auto Echo Reducer	•	-	-	-	-
Reverb Reducer	•	•	-	-	-
<b>Sound Quality &amp; Control</b>					
Bandwidth extension	•	-	-	-	-
Bandwidth compression	•	•	•	•	•
Feedback preventer	•	•	•	•	•
Music Enhancer (presets)	3	3	1	1	-
Automatic Acclimatization	•	•	•	•	•
<b>Audio Streaming</b>					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
ASHA	-	-	-	-	-
<b>Headset Mode</b>					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
<b>TeleCare</b>					
Full live remote tuning	-	-	-	-	-
Basic portal remote tuning	•	•	•	•	•
<b>Accessories &amp; Apps</b>		<b>Chargers</b>		<b>CROS</b>	
Smart Key	•	Charging Station	-	Reach CROS R-Li T BC	-
TV Sound	-	Travel Charger	-	Reach CROS R-Li	-
Smart Transmitter 2,4	-	Charging+ Station	-	Reach CROS SR	-
Smart Mic	-	IF CIC Travel Charger	-	Reach CROS iX-CIC Li	•
Rexton App	•	Travel Charger Custom CIC	•		
<b>Programming</b>		Universal Cable Charger			
Noahlink Wireless	-				
ConnexxAir	•				
HI-PRO	-				