

REACH Custom IIC AND CIC

A NEXT-LEVEL CUSTOM MADE INVISIBLE

The Reach Custom IIC and CIC are designed to optimize the perception of small size in several ways and with our new I3 faceplate the invisibility of these devices are next level.




- > Our smallest customized solution
- > 36% better battery lifetime compared to predecessor (~75h)
- > Ultra HD e2e binaural link enables OneMic Directionality
- > Rexton App compatible
- > Programmable via ConnexxAir and/or programming adaptor
- > Lifeproof

Colors



Fitting range is based on 2 ccm coupler:



* Available for both IIC and CIC

** Available for CIC

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	36/18/6	32/16/6	24/12/4	16/8/4
Battery	10	10	10	10	10
T-coil	-	-	-	-	-
Bluetooth	-	-	-	-	-
Tinnitus function (Sound therapy and Notch)	•	•	•	•	-
Wireless sync	•	•	•	•	•
IP68	-	-	-	-	-
ATEX	-	-	-	-	-
EarWear	3.0	3.0	3.0	3.0	3.0
Technology Highlights					
Soundpro 3.0	•	•	•	•	•
Multi-Voice Focus	-	-	-	-	-
Focus Streams	-	-	-	-	-
Voice Stabilizer	•	•	•	•	•
Speech Preservation Technology	-	-	-	-	-
My Voice 2.0	-	-	-	-	-
Motion Sensor for scene analysis	-	-	-	-	-
Directionality					
OneMic directionality	•	•	•	-	-
Automatic adaptive	•	•	•	•	•
iOmni Sound Locator	-	-	-	-	-
iFocus 360 (Front/Back)	-	-	-	-	-
iFocus 360 (Left /Right)	-	-	-	-	-
Auto Stereo iLock (Narrow)	-	-	-	-	-
Noise Reduction					
Noise Management	•	•	•	•	•
Sound Smoothing	•	•	•	•	•
Wind Noise Reduction	•	•	•	•	•
Auto Echo Reducer	•	-	-	-	-
Reverb Reducer	Manual	Manual	-	-	-
Sound Quality & Control					
Bandwidth extension	•	-	-	-	-
Bandwidth compression	•	•	•	•	•
Feedback preventer	•	•	•	•	•
Music Enhancer (presets)	3	3	1	1	-
Smart Automatic Equalizer	•	•	•	•	•
Smart Automatic Acclimatization	•	•	•	•	•
Automatic Classifier	•	•	•	•	•
Data Logging	•	•	•	•	•
Audio Streaming					
BT LE Audio	-	-	-	-	-
Made for iPhone	-	-	-	-	-
ASHA (Android 10 or higher)	-	-	-	-	-
Headset Mode					
iOS	-	-	-	-	-
BT LE Audio	-	-	-	-	-
App					
Rexton App	•	•	•	•	•
Spatial configurator	-	-	-	-	-
Rexton Assist	-	-	-	-	-
TeleCare					
Full live remote tuning	-	-	-	-	-
Basic portal remote tuning	•	•	•	•	•
Accessories			CROS		
Smart Key	•		Reach CROS R-Li	-	
Smart Transmitter 2,4	-		Reach CROS SR	-	
Smart Mic	-		Reach CROS iX-CIC Li	•	
Fitting					
Noahlink Wireless	-				
ConnexxAir	•				
Programming adaptor	•				

• available - not available