

REACH CUSTOM ITC AND ITE




- › Enabling users to be fully involved in group conversations with Multi-Voice Focus
- › With Voice Stabilizer the contrast in every conversation is controlled for easier speech intelligibility
- › Battery driven solutions (size 312) for a full days use:
 - › Custom ITC and ITE delivers up to 135h runtime depending on receiver strength and microphone set-up
- › Lifeproof to withstand everyday challenges


Faceplate configuration:

- › N faceplate:
 - › One-mic configuration available for ITC
- › N4 faceplate:
 - › Two-mic configuration available for ITC and ITE
- › Both faceplates offers:
 - › Telecoil
 - › Push button and/or volume control wheel
 - › Non-wireless option

Colors



Fitting range is based on 2 ccm coupler:



- | | | |
|---|---|---|
| • 50 dB / 113 dB SPL
(2 ccm coupler) | • 55 dB / 118 dB SPL
(2 ccm coupler) | • 60 dB / 124 dB SPL
(2 ccm coupler) |
| • 60 dB / 124 dB SPL
(Ear simulator) | • 65 dB / 127 dB SPL
(Ear simulator) | • 70 dB / 133 dB SPL
(Ear simulator) |

* Available for ITC

** Available for both ITC and ITE

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	36/18/6	32/16/6	24/12/4	16/8/4
LED	-	-	-	-	-
T-coil	•	•	•	•	•
Bluetooth	-	-	-	-	-
Tinnitus Sound function	•	•	•	•	-
Tinnitus Notch function	•	•	•	•	-
Wireless sync	•	•	•	•	•
IP68	-	-	-	-	-
XPhone	•	•	•	•	•
Technology Highlights					
Multi-Voice Focus	•	•	•	-	-
Focus Beams	4	2	1	1	1
Voice Stabilizer	•	•	•	•	•
Speech Preservation Technology	•	•	•	•	•
My Voice 2.0	-	-	-	-	-
Motion sensor for scene analysis	-	-	-	-	-
Directionality					
Automatic adaptive	•	•	•	•	•
iOmni Sound Locator	•	•	•	•	•
Auto iFocus 360	•	-	-	-	-
iFocus 360	•	•	-	-	-
Auto Focus 360	-	•	-	-	-
Auto Stereo iLock	•	•	•	-	-
OneMic Directionality	-	-	-	-	-
Noise Reduction					
Noise Management	•	•	•	•	•
Sound Smoothing	•	•	•	•	•
Wind Noise Reduction	•	•	•	•	•
Auto Echo Reducer	•	-	-	-	-
Reverb Reducer	•	•	-	-	-
Sound Quality & Control					
Bandwidth extension	•	-	-	-	-
Bandwidth compression	•	•	•	•	•
Feedback preventer	•	•	•	•	•
Music Enhancer (presets)	3	3	1	1	-
Automatic Acclimatization	•	•	•	•	•
Audio Streaming					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
ASHA	-	-	-	-	-
Headset Mode					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
TeleCare					
Full live remote tuning	-	-	-	-	-
Basic portal remote tuning	•	•	•	•	•
Accessories & Apps			CROS		
Smart Key	•		Reach CROS R-Li T BC	-	
TV Sound	-		Reach CROS R-Li	-	
Smart Transmitter 2.4	-		Reach CROS SR	-	
Smart Mic	-		Reach CROS iX-CIC Li	•	
Rexton App	•				
Programming					
Noahlink Wireless	-				
ConnexxAir	•				
HI-PRO	•				

REACH CUSTOM IIC AND CIC

A NEXT-LEVEL CUSTOM MADE INVISIBLE



- › With Voice Stabilizer the contrast in every conversation is controlled for easier speech intelligibility
- › OneMic Directionality enables the benefits of a directional microphone
- › Our smallest customized solution yet
- › Battery driven solutions with 36% better battery lifetime compared to predecessor:
 - › Custom IIC and CIC (battery size 10) delivers up to 75h runtime
- › Lifeproof to withstand everyday challenges

Colors

Black




Brown



Mocha



Fitting range is based on 2 ccm coupler:



* Available for both IIC and CIC

** Available for CIC

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	36/18/6	32/16/6	24/12/4	16/8/4
LED	-	-	-	-	-
T-coil	-	-	-	-	-
Bluetooth	-	-	-	-	-
Tinnitus Sound function	•	•	•	•	-
Tinnitus Notch function	•	•	•	•	-
Wireless sync	•	•	•	•	•
IP68	-	-	-	-	-
XPhone	-	-	-	-	-
Technology Highlights					
Multi-Voice Focus	-	-	-	-	-
Focus Beams	-	-	-	-	-
Voice Stabilizer	•	•	•	•	•
Speech Preservation Technology	-	-	-	-	-
My Voice 2.0	-	-	-	-	-
Motion sensor for scene analysis	-	-	-	-	-
Directionality					
Automatic adaptive	•	•	•	•	•
iOmni Sound Locator	-	-	-	-	-
Auto iFocus 360	-	-	-	-	-
iFocus 360	-	-	-	-	-
Auto Focus 360	-	-	-	-	-
Auto Stereo iLock	-	-	-	-	-
OneMic Directionality	•	•	•	-	-
Noise Reduction					
Noise Management	•	•	•	•	•
Sound Smoothing	•	•	•	•	•
Wind Noise Reduction	•	•	•	•	•
Auto Echo Reducer	•	-	-	-	-
Reverb Reducer	•	•	-	-	-
Sound Quality & Control					
Bandwidth extension	•	-	-	-	-
Bandwidth compression	•	•	•	•	•
Feedback preventer	•	•	•	•	•
Music Enhancer (presets)	3	3	1	1	-
Automatic Acclimatization	•	•	•	•	•
Audio Streaming					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
ASHA	-	-	-	-	-
Headset Mode					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
TeleCare					
Full live remote tuning	-	-	-	-	-
Basic portal remote tuning	•	•	•	•	•
Accessories & Apps			CROS		
Smart Key	•		Reach CROS R-Li T BC	-	
TV Sound	-		Reach CROS R-Li	-	
Smart Transmitter 2.4	-		Reach CROS SR	-	
Smart Mic	-		Reach CROS iX-CIC Li	•	
Rexton App	•				
Programming					
Noahlink Wireless	-				
ConnexxAir	•				
HI-PRO	•				