


REACH INOX-CIC LI

THE FIRST INVISIBLE HEARING AID WITH RECHARGEABILITY



- > With Voice Stabilizer the contrast in every conversation is controlled for easier speech intelligibility
- > OneMic Directionality enables the benefits of a directional microphone
- > A discreet design with easy handling makes it virtually invisible
- > Portable charger with built-in power bank
- > Rechargeable solution for a full days use:
 - > inoX-CIC Li delivers up to 28h runtime
- > Lifeproof to withstand everyday challenges
- > ATEX certified

Faceplate colors



- 50 dB / 114 dB SPL (2 ccm coupler)
- 61 dB / 126 dB SPL (Ear simulator)



For a customized, yet invisible rechargeable fit, Reach inoX-CIC Li is compatible to earmold fittings.

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	36/18/6	32/16/6	24/12/4	16/8/4
LED	-	-	-	-	-
T-coil	-	-	-	-	-
Bluetooth	-	-	-	-	-
Tinnitus Sound function	•	•	•	•	-
Tinnitus Notch function	•	•	•	•	-
Wireless sync	•	•	•	•	•
IP68	•	•	•	•	•
XPhone	•	•	•	•	•
Technology Highlights					
Multi-Voice Focus	-	-	-	-	-
Focus Beams	-	-	-	-	-
Voice Stabilizer	•	•	•	•	•
Speech Preservation Technology	-	-	-	-	-
My Voice 2.0	-	-	-	-	-
Motion sensor for scene analysis	-	-	-	-	-
Directionality					
Auto Dir	•	•	•	-	-
Directional Adaptive	•	•	•	-	-
iOmni Sound Locator	-	-	-	-	-
Auto iFocus 360	-	-	-	-	-
iFocus 360	-	-	-	-	-
Auto Focus 360	-	-	-	-	-
Auto Stereo iLock	-	-	-	-	-
OneMic Directionality	•	•	•	-	-
Noise Reduction					
Noise Management	•	•	•	•	•
Sound Smoothing	•	•	•	•	•
Wind Noise Reduction	•	•	•	•	•
Auto Echo Reducer	•	-	-	-	-
Reverb Reducer	•	•	-	-	-
Sound Quality & Control					
Bandwidth extension	•	-	-	-	-
Bandwidth compression	•	•	•	•	•
Feedback preventer	•	•	•	•	•
Music Enhancer (presets)	3	3	1	1	-
Automatic Acclimatization	•	•	•	•	•
Audio Streaming					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
ASHA	-	-	-	-	-
Headset Mode					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
TeleCare					
Full live remote tuning	-	-	-	-	-
Basic portal remote tuning	•	•	•	•	•
Accessories & Apps		Chargers		CROS	
Remote Control	•	Charging Station	-	Reach CROS R-Li T BC	-
TV Sound	-	Travel Charger	•	Reach CROS R-Li	-
Smart Mic	-	Charging+ Station	-	Reach CROS SR	-
Rexton App	•	Universal Cable Charger	-	Reach CROS iX-CIC Li	•
Programming					
Noahlink Wireless	-				
ConnexxAir	•				
HI-PRO	-				