

RECHARGE TECHNOLOGY



REXTON ReCharge will bring new energy to hearing situations everywhere. Bringing proven hearing aid technology and will feature a handy cable charger for recharging at home or on the move. No more hassles with tiny batteries. Just plug in, power up and away you go.

EASY TO USE

No more hassle with tiny batteries since ReCharge is rechargeable. They are easy to fit and programmable with Noahlink wireless design.

AFFORDABILITY

The first rechargeable hearing aid with proven hearing technology in the essential portfolio. The charger is designed for optimal usability and affordability.

RECHARGEABILITY

With a unique performance level offering, you can recharge at home or on the go. Plug in and power up with the handy cable-charger which can recharge one or two hearing aids.

LIFEPROOF

REXTON hearing aids are built using robust components and tested extensively to ensure they can withstand whatever life throws at them. They are durable enough to withstand everyday challenges and harsh environments. All ReCharge products are IP68 certified and ATEX.

TECHNOLOGY

- › Localization in quieter environments is improved by mimicking the ear acoustics.
- › Extended Dynamic Range
- › High input sounds with improved sound quality for less distortion.
- › Adaptive compression for more natural speech.
- › 48-channel first fit for a more precise and easy fitting from scratch



RECHARGE R-LI

JUST PLUG IN AND POWER UP

The simplicity and affordability of REXTON ReCharge will have a major impact. New energy for work, hobbies, travel, nights out, for time with the family, and for conversations everywhere. REXTON ReCharge is powered by proven technology. No more hassles with handling and changing tiny batteries. Just plug in and power up, and everything's ready to go.




- > Programmable via Noahlink Wireless
- > Takes you through the day and night with 35H runtime
- > Unique rechargeability offering in the 30,20,10,5 performance levels
- > Addresses mild to severe hearing loss
- > Water-resistant, IP68 rating
- > Lifeproof
- > ATEX

Colors



All fitting range is based on 2 ccm coupler:



Features / Performance Level	30	20	10	5
Channels/Controls/Programs	16/8/4	16/8/4	12/6/4	8/4/4
Battery	Lithium-ion	Lithium-ion	Lithium-ion	Lithium-ion
Tcoil	-	-	-	-
Tinnitus Sound therapy	•	-	-	-
Wireless sync	•	•	-	-
IP68	•	•	•	•
ATEX	•	•	•	•
EarWear	miniReceiver 2.0 + click dome/sleeve/Earmold			
Technology highlights				
First Fit 48 channels	•	•	•	•
Voice Ranger	•	•	-	-
Adaptive compression	•	•	•	•
Extended Dynamic Range	•	•	•	•
Directionality				
Multichannel Adaptive Directional Microphone	•	-	-	-
Automatic Directional Microphone	•	•	•	-
Fixed Directional Microphone	•	•	•	•
Omni Sound Locator	•	•	-	-
Noise Reduction				
Noise Management	•••	••	••	•
Sound Smoothing	•••	-	-	-
Wind Noise Cancellation	•	-	-	-
Sound quality & control				
Bandwidth compression	•	•	-	-
Feedback preventer	•••	••	••	••
Smart Automatic Equalizer	•	•	-	-
Smart Automatic Acclimatization	•••	••	-	-
Automatic Classifier	•••	•	•	•
Data Logging	•••	•••	••	••
App				
Rexton App	•	•	•	•
TeleCare				
Full live remote tuning	-	-	-	-
Basic remote tuning	•	•	•	•

Accessories	
Smart Key	o
Smart Transmitter 2,4	-
Smart Mic	-
Fitting	
Noahlink Wireless	mandatory
Connexx Air	-

Chargers	
Charging Station R	-
Travel Charger RIC / Slim RIC	-
Charging+ Station (with Dry&Clean function)	-
Instant Fit CIC Travel Charger	-
Cable Charger	mandatory



RECHARGE B-LI P


JUST PLUG IN AND POWER UP

The simplicity and affordability of REXTON ReCharge B-Li P will secure new energy for an entire day. Moderate to severe hearing losses can benefit from proven technology for clear and comfortable sounds.




- › Programmable via Noahlink Wireless
- › Takes you through the day and night with 39H runtime
- › Unique rechargeability offering in the 30,20,10,5 performance levels
- › Addresses mild to severe hearing loss
- › Water-resistant, IP68 rating
- › Lifeproof
- › ATEX

Colors



All fitting range is based on 2 ccm coupler:



- 77 dB / 135 dB SPL (2 ccm coupler)

Features / Performance Level	30	20	10	5
Channels/Controls/Programs	16/8/4	16/8/4	12/6/4	8/4/4
Battery	Lithium-ion	Lithium-ion	Lithium-ion	Lithium-ion
Tcoil	ReCharge B-Li P T	ReCharge B-Li P T	-	-
Tinnitus Sound therapy	•	-	-	-
Wireless sync	•	•	-	-
IP68	•	•	•	•
ATEX	•	•	•	•
EarWear	Earmold only (no thin tube available)			
Technology highlights				
First Fit 48 channels	•	•	•	•
Voice Ranger	•	•	-	-
Adaptive compression	•	•	•	•
Extended Dynamic Range	•	•	•	•
Directionality				
Multichannel Adaptive Directional Microphone	•	-	-	-
Automatic Directional Microphone	•	•	•	-
Fixed Directional Microphone	•	•	•	•
Omni Sound Locator	•	•	-	-
Noise Reduction				
Noise Management	•••	••	••	•
Sound Smoothing	•••	-	-	-
Wind Noise Reduction	•	-	-	-
Sound quality & control				
Bandwidth compression	•	•	-	-
Feedback preventer	•••	••	••	••
Smart Automatic Equalizer	•	•	-	-
Smart Automatic Acclimatization	•••	••	-	-
Automatic Classifier	••••	•	•	•
Data Logging	••••	••••	••	••
App				
Rexton App	•	•	•	•
TeleCare				
Full live remote tuning	-	-	-	-
Basic remote tuning	•	•	•	•

Accessories	
Smart Key	o
Smart Transmitter 2,4	-
Smart Mic	-
Fitting	
Noahlink Wireless	mandatory
Connexx Air	-

Chargers	
Charging Station R	-
Travel Charger RIC / Slim RIC	-
Charging+ Station (with Dry&Clean function)	-
Instant Fit CIC Travel Charger	-
Cable Charger	mandatory



RECHARGE B-LI HP


JUST PLUG IN AND POWER UP

The powerful long-lasting REXTON ReCharge B-Li HP will secure new energy for up to 72h. Severe to profound hearing losses can benefit from proven technology for clear and comfortable sounds.




- > Takes you through the day and night with 88H runtime
- > Unique rechargeability offering in the 200, 100, 75 performance levels
- > Addresses severe to profound hearing loss
- > Programmable via Noahlink Wireless
- > Water-resistant, IP68 rating
- > Lifeproof
- > ATEX

Colors



All fitting range is based on 2 ccm coupler:



- 82 dB / 140 dB SPL (2 ccm coupler)

Features / Performance Level	30	20	10	5
Channels/Controls/Programs	16/8/4	16/8/4	12/6/4	8/4/4
Battery	Lithium-ion	Lithium-ion	Lithium-ion	Lithium-ion
Tcoil	-	-	-	-
Tinnitus Sound therapy	•	-	-	-
Wireless sync	•	•	-	-
IP68	•	•	•	•
ATEX	•	•	•	•
EarWear	Earmold only (no thin tube available)			

Technology highlights				
First Fit 48 channels	•	•	•	•
Voice Ranger	•	•	-	-
Adaptive compression	•	•	•	•
Extended Dynamic Range	•	•	•	•

Directionality				
Multichannel Adaptive Directional Microphone	•	-	-	-
Automatic Directional Microphone	•	•	•	-
Fixed Directional Microphone	•	•	•	•
Omni Sound Locator	•	•	-	-

Noise Reduction				
Noise Management	•••	••	••	•
Sound Smoothing	•••	-	-	-
Wind Noise Reduction	•	-	-	-

Sound quality & control				
Bandwidth compression	•	•	-	-
Feedback preventer	•••	••	••	••
Smart Automatic Equalizer	•	•	-	-
Smart Automatic Acclimatization	•••	••	-	-
Automatic Classiher	••••	•	•	•
Data Logging	••••	••••	••	••

App				
Rexton App	•	•	•	•

TeleCare				
Full live remote tuning	-	-	-	-
Basic remote tuning	•	•	•	•

Accessories	
Smart Key	o
Smart Transmitter 2,4	-
Smart Mic	-

Chargers	
Charging Station R	-
Travel Charger RIC / Slim RIC	-
Charging+ Station (with Dry&Clean function)	-
Instant Fit CIC Travel Charger	-
Cable Charger	mandatory

Fitting	
Noahlink Wireless	mandatory
Connexx Air	-