REACH B-LI RUGGED

OUR TOUGHEST HEARING AID NOW WITH MULTI-VOICE FOCUS




- Enabling users to be fully involved in group conversations with Multi-Voice Focus
- With My Voice 2.0 the users own voice stays comfortable with a dedicated processing separated from the target speech sounds
- Convenient window in the housing for display of serial number
- Hands-free phone calls via Headset Mode with Call Control when connected via MFi or LE Audio
- Rechargeable solution for a full days use:
  - Rugged delivers up to 32h runtime (incl. 5h streaming)
  - Rugged delivers up to 36h runtime without streaming
- Lifeproof to withstand everyday challenges
- AI with Rexton Assist helps you to seamlessly customize your hearing aids

Rugged offers a uniquely tough build:

- Water resistant construction
- Sweat proof
- Drop proof
- Scratch resistant
- IP68 certified
- ATEX certified

Colors





Earhook

- 60 dB / 133 dB SPL (2 ccm coupler)
- 67 dB / 138 dB SPL (Ear simulator)

TOUGH DESIGN WITH A SUPER-COMPACT SIZE

Our toughest hearing aid is now available within our platform Reach. Reach Rugged will deliver excellent sound processing in a small and secure package. It's structure has an airtight receiver, sealed charging contacts, reinforced receiver suspension, thickened housing, and sealed parting lines.

Size comparison



Features / Performance Level	80	60	40
Channels/Controls/Programs	48/20/6	36/18/6	32/16/6
LED	-	-	-
T-coil	•	•	•
Bluetooth	•	•	•
Tinnitus Sound function	•	•	•
Tinnitus Notch function	•	•	•
Wireless sync	•	•	•
IP68	•	•	•
XPhone	•	•	•
Technology Highlights			
Multi-Voice Focus	•	•	•
Focus Beams	4	2	1
Voice Stabilizer	•	•	•
Speech Preservation Technology	•	•	•
My Voice 2.0	•	•	•
Motion sensor for scene analyzis	•	•	•
Directionality			
Auto Dir	•	•	•
Directional Adaptive	•	•	•
iOmni Sound Locator	•	•	•
Auto iFocus 360	•	-	-
iFocus 360	•	•	-
Auto Focus 360	-	•	-
Auto Stereo iLock	•	•	•
OneMic Directionality	-	-	-
Noise Reduction			
Noise Management	•	•	•
Sound Smoothing	•	•	•
Wind Noise Reduction	•	•	•
Auto Echo Reducer	•	-	-
Reverb Reducer	•	•	-
Sound Quality & Control			
Bandwidth extension	•	-	-
Bandwidth compression	•	•	•
Feedback Preventer	•	•	•
Music Enhancer (presets)	3	3	1
Automatic Acclimatization	•	•	•
Audio Streaming			
Bluetooth Classic	-	-	-
LE Audio	•	•	•
MFi	•	•	•
ASHA	•	•	•
Headset Mode			
Bluetooth Classic	-	-	-
LE Audio	•	•	•
MFi	•	•	•
TeleCare			
Full live remote tuning	•	•	•
Basic portal remote tuning	•	•	•

Accessories & Apps	
Remote Control	•
TV Sound	•
Smart Mic	•
Rexton App	•

Programming	
Noahlink Wireless	•
ConnexxAir	-
HI-PRO	-

Chargers	
Charging Station	•
Travel Charger	-
Charging+ Station	-
Universal Cable Charger	•

CROS	
Reach CROS R-Li T BC	-
Reach CROS R-Li	•
Reach CROS SR	-
Reach CROS iX-CIC Li	-

• available - not available