
REACH B-LI P AND B-LI HP

› Enabling users to be fully involved in group conversations with Multi-Voice Focus
› With My Voice 2.0, the sound of the user's own voice stays comfortable with a dedicated processing separated from the target speech sounds
› Convenient window in the housing for display of serial number
› Hands-free phone calls via Headset Mode with Call Control when connected via MFi or LE Audio
› Tamper-proof rocker switch and earhook
› Lifeproof to withstand everyday challenges
› ATEX certified*
› AI with Rexton Assist helps you to seamlessly customize your hearing aids


Colors



REACH B-LI P

› Rechargeable solution for a full days use:
  › B-Li P delivers up to 30h runtime (incl. 5h of streaming) or up to 35h (without streaming).






Earhook

- 77 dB / 135 dB SPL (2 ccm coupler)
- 82 dB / 139 dB SPL (Ear simulator)

REACH B-LI HP

› Rechargeable solution for a full days use:
  › B-Li HP delivers up to 65h runtime (incl. 5h streaming) or up to 73h (without streaming).





Earhook

- 82 dB / 140 dB SPL (2 ccm coupler)
- 85 dB / 143 dB SPL (Ear simulator)

* The ATEX certificate is planned to be available within six months after launch.

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6[1]	36/18/6	32/16/6	24/12/4	16/8/4
LED	•	•	•	-	-
T-coil	•	•	•	•	•
Bluetooth	•	•	•	•	•
Tinnitus Sound function	•	•	•	•	-
Tinnitus Notch function[2]	•	•	•	•	-
Wireless sync	•	•	•	•	•
IP68	•	•	•	•	•
XPhone	•	•	•	•	•
Technology Highlights					
Multi-Voice Focus	•	•	•	-	-
Focus Beams	4	2	1	-	-
Voice Stabilizer	•	•	•	•	•
Speech Preservation Technology	•	•	•	•	•
My Voice 2.0	•	•	•	-	-
Motion sensor for scene analyzis	•	•	•	-	-
Directionality					
Auto Dir	•	•	•	•	•
Directional Adaptive	•	•	•	•	-
iOmni Sound Locator	•	•	•	•	•
Auto iFocus 360	•	-	-	-	-
iFocus 360	•	•	-	-	-
Auto Focus 360	-	•	-	-	-
Auto Stereo iLock	•	•	•	-	-
OneMic Directionality	-	-	-	-	-
Noise Reduction					
Noise Management	•	•	•	•	•
Sound Smoothing	•	•	•	•	•
Wind Noise Reduction	•	•	•	•	•
Auto Echo Reducer	•	-	-	-	-
Reverb Reducer	•	•	-	-	-
Sound Quality & Control					
Bandwidth extension[2]	•	-	-	-	-
Bandwidth compression	•	•	•	•	•
Feedback preventer	•	•	•	•	•
Music Enhancer (presets)	3	3	1	1	-
Automatic Acclimatization	•	•	•	•	•
Audio Streaming					
Bluetooth Classic	-	-	-	-	-
LE Audio	•	•	•	•	•
MFi	•	•	•	•	•
ASHA	•	•	•	•	•
Headset Mode					
Bluetooth Classic	-	-	-	-	-
LE Audio	•	•	•	•	•
MFi	•	•	•	•	•
TeleCare					
Full live remote tuning	•	•	•	•	•
Basic portal remote tuning	•	•	•	•	•

Accessories & Apps	
Remote Control	•
TV Sound	•
Smart Mic	•
Rexton App	•

Programming	
Noahlink Wireless	•
ConnexxAir	-
HI-PRO	-

Chargers	
Charging Station	•
Travel Charger	-
Charging+ Station	•
Universal Cable Charger	•

CROS	
Reach CROS R-Li T BC	-
Reach CROS R-Li	•
Reach CROS SR	-
Reach CROS iX-CIC Li	-

1) 34/16/6 for B-Li HP
2) Not available for B-Li HP

• available - not available