REACH CUSTOM ITC AND ITE



› Enabling users to be fully involved in group conversations with Multi-Voice Focus

› With Voice Stabilizer the contrast in every conversation is controlled for easier speech intelligibility

› Battery driven solutions (size 312) for a full days use:

  › Custom ITC and ITE delivers up to 135h runtime depending on receiver strength and microphone set-up

› Lifeproof to withstand everyday challenges

Faceplate configuration:


› N faceplate:
  › One-mic configuration available for ITC
› N4 faceplate:
  › Two-mic configuration available for ITC and ITE

› Both faceplates offers optional:
  › Telecoil
  › Push button and/or volume control wheel
  › Non-wireless option

Colors


Black Tan Beige Mocha

S-Receiver*




- 50 dB / 113 dB SPL (2 ccm coupler)
- 60 dB / 124 dB SPL (Ear simulator)

M-Receiver**



- 55 dB / 118 dB SPL (2 ccm coupler)
- 65 dB / 127 dB SPL (Ear simulator)

P-Receiver**



- 60 dB / 124 dB SPL (2 ccm coupler)
- 70 dB / 133 dB SPL (Ear simulator)

* Available for ITC
** Available for both ITC and ITE

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	36/18/6	32/16/6	24/12/4	16/8/4
LED	-	-	-	-	-
T-coil	o	o	o	o	o
Bluetooth	-	-	-	-	-
Tinnitus Sound function	•	•	•	•	-
Tinnitus Notch function	•	•	•	•	-
Wireless sync	•	•	•	•	•
IP68	-	-	-	-	-
XPhone	•	•	•	•	•
Technology Highlights					
Multi-Voice Focus	•	•	•	-	-
Focus Beams	4	2	1	1	1
Voice Stabilizer	•	•	•	•	•
Speech Preservation Technology	•	•	•	•	•
My Voice 2.0	-	-	-	-	-
Motion sensor for scene analyzis	-	-	-	-	-
Directionality					
Auto Dir	•	•	•	•	•
Directional Adaptive	•	•	•	•	-
iOmni Sound Locator	-	-	-	-	-
Auto iFocus 360	•	-	-	-	-
iFocus 360	•	•	-	-	-
Auto Focus 360	-	•	-	-	-
Auto Stereo iLock	•	•	•	-	-
OneMic Directionality	-	-	-	-	-
Noise Reduction					
Noise Management	•	•	•	•	•
Sound Smoothing	•	•	•	•	•
Wind Noise Reduction	•	•	•	•	•
Auto Echo Reducer	•	-	-	-	-
Reverb Reducer	•	•	-	-	-
Sound Quality & Control					
Bandwidth extension	•	-	-	-	-
Bandwidth compression	•	•	•	•	•
Feedback preventer	•	•	•	•	•
Music Enhancer (presets)	3	3	1	1	-
Automatic Acclimatization	•	•	•	•	•
Audio Streaming					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
ASHA	-	-	-	-	-
Headset Mode					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
TeleCare					
Full live remote tuning	-	-	-	-	-
Basic portal remote tuning	•	•	•	•	•

Accessories & Apps	
Remote Control	•
TV Sound	-
Smart Mic	-
Rexton App	•

CROS	
Reach CROS R-Li T BC	-
Reach CROS R-Li	-
Reach CROS SR	-
Reach CROS iX-CIC Li	•

Programming	
Noahlink Wireless	-
ConnexxAir	•
HI-PRO	•

Features are based on a two-mic configuration with wireless option.

• available - not available o optional

REACH CUSTOM IIC AND CIC

A NEXT-LEVEL CUSTOM MADE INVISIBLE




- With Voice Stabilizer the contrast in every conversation is controlled for easier speech intelligibility
- OneMic Directionality enables the benefits of a directional microphone
- Our smallest customized solution yet
- Battery driven solutions with 36% better battery lifetime compared to predecessor:
  - Custom IIC and CIC (battery size 10) delivers up to 75h runtime
- Lifeproof to withstand everyday challenges

Colors

Black


Mocha

S-Receiver*




- 50 dB / 113 dB SPL (2 ccm coupler)
- 60 dB / 124 dB SPL (Ear simulator)

M-Receiver**



- 55 dB / 118 dB SPL (2 ccm coupler)
- 65 dB / 127 dB SPL (Ear simulator)

P-Receiver**



- 60 dB / 124 dB SPL (2 ccm coupler)
- 70 dB / 133 dB SPL (Ear simulator)

* Available for both IIC and CIC
** Available for CIC

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	36/18/6	32/16/6	24/12/4	16/8/4
LED	-	-	-	-	-
T-coil	-	-	-	-	-
Bluetooth	-	-	-	-	-
Tinnitus Sound function	•	•	•	•	-
Tinnitus Notch function	•	•	•	•	-
Wireless sync	•	•	•	•	•
IP68	-	-	-	-	-
XPhone	•	•	•	•	•
Technology Highlights					
Multi-Voice Focus	-	-	-	-	-
Focus Beams	-	-	-	-	-
Voice Stabilizer	•	•	•	•	•
Speech Preservation Technology	-	-	-	-	-
My Voice 2.0	-	-	-	-	-
Motion sensor for scene analyzis	-	-	-	-	-
Directionality					
Auto Dir	•	•	•	-	-
Directional Adaptive	•	•	•	-	-
iOmni Sound Locator	-	-	-	-	-
Auto iFocus 360	-	-	-	-	-
iFocus 360	-	-	-	-	-
Auto Focus 360	-	-	-	-	-
Auto Stereo iLock	-	-	-	-	-
OneMic Directionality	•	•	•	-	-
Noise Reduction					
Noise Management	•	•	•	•	•
Sound Smoothing	•	•	•	•	•
Wind Noise Reduction	•	•	•	•	•
Auto Echo Reducer	•	-	-	-	-
Reverb Reducer	•	•	-	-	-
Sound Quality & Control					
Bandwidth extension	•	-	-	-	-
Bandwidth compression	•	•	•	•	•
Feedback preventer	•	•	•	•	•
Music Enhancer (presets)	3	3	1	1	-
Automatic Acclimatization	•	•	•	•	•
Audio Streaming					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
ASHA	-	-	-	-	-
Headset Mode					
Bluetooth Classic	-	-	-	-	-
LE Audio	-	-	-	-	-
MFi	-	-	-	-	-
TeleCare					
Full live remote tuning	-	-	-	-	-
Basic portal remote tuning	•	•	•	•	•

Accessories & Apps	
Remote Control	•
TV Sound	-
Smart Mic	-
Rexton App	•

CROS	
Reach CROS R-Li T BC	-
Reach CROS R-Li	-
Reach CROS SR	-
Reach CROS iX-CIC Li	•

Programming	
Noahlink Wireless	-
ConnexxAir	•
HI-PRO	•

Features are based on wireless option.

• available - not available