
REACH R-LI AND R-LI T




- › Enabling users to be fully involved in group conversations with Multi-Voice Focus
- › Hands-free phone calls via Headset Mode with Call Control when connected via MFi or LE Audio
- › Portable charger with built-in power bank
- › Rechargeable solutions for a full days use:
  - › R-Li delivers up to 24h runtime (incl. 5h streaming)
  - › R-Li T delivers up to 34h runtime (incl. 5h streaming)
- › Lifeproof to withstand everyday challenges
- › ATEX certified
- › AI with Rexton Assist helps you to seamlessly customize your hearing aids

Colors




S-Receiver




- 46 dB / 110 dB SPL (2 ccm coupler)
- 56 dB / 120 dB SPL (Ear simulator)

M-Receiver




- 60 dB / 119 dB SPL (2 ccm coupler)
- 70 dB / 129 dB SPL (Ear simulator)

P-Receiver



- 65 dB / 122 dB SPL (2 ccm coupler)
- 75 dB / 131 dB SPL (Ear simulator)

HP-Receiver



- 75 dB / 131 dB SPL (2 ccm coupler)
- 83 dB / 138 dB SPL (Ear simulator)

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	36/18/6	32/16/6	24/12/4	16/8/4
LED	-	-	-	-	-
T-coil	R-Li T	R-Li T	R-Li T	R-Li T	R-Li T
Bluetooth	•	•	•	•	•
Tinnitus Sound function	•	•	•	•	-
Tinnitus Notch function	•	•	•	•	-
Wireless sync	•	•	•	•	•
IP68	•	•	•	•	•
XPhone	•	•	•	•	•
Technology Highlights					
Multi-Voice Focus	•	•	•	-	-
Focus Beams	4	2	1	-	-
Voice Stabilizer	•	•	•	•	•
Speech Preservation Technology	•	•	•	•	•
My Voice 2.0	•	•	•	-	-
Motion sensor for scene analyzis	•	•	•	-	-
Directionality					
Auto Dir	•	•	•	•	•
Directional Adaptive	•	•	•	•	-
iOmni Sound Locator	•	•	•	•	•
Auto iFocus 360	•	-	-	-	-
iFocus 360	•	•	-	-	-
Auto Focus 360	-	•	-	-	-
Auto Stereo iLock	•	•	•	-	-
OneMic Directionality	-	-	-	-	-
Noise Reduction					
Noise Management	•	•	•	•	•
Sound Smoothing	•	•	•	•	•
Wind Noise Reduction	•	•	•	•	•
Auto Echo Reducer	•	-	-	-	-
Reverb Reducer	•	•	-	-	-
Sound Quality & Control					
Bandwidth extension	•	-	-	-	-
Bandwidth compression	•	•	•	•	•
Feedback preventer	•	•	•	•	•
Music Enhancer (presets)	3	3	1	1	-
Automatic Acclimatization	•	•	•	•	•
Audio Streaming					
Bluetooth Classic	-	-	-	-	-
LE Audio	•	•	•	•	•
MFi	•	•	•	•	•
ASHA	•	•	•	•	•
Headset Mode					
Bluetooth Classic	-	-	-	-	-
LE Audio	•	•	•	•	•
MFi	•	•	•	•	•
TeleCare					
Full live remote tuning	•	•	•	•	•
Basic portal remote tuning	•	•	•	•	•

Accessories & Apps	
Remote Control	•
TV Sound	•
Smart Mic	•
Rexton App	•

Chargers	
Charging Station	•
Travel Charger	•
Charging+ Station	•
Universal Cable Charger	•

CROS	
Reach CROS R-Li T BC	-
Reach CROS R-Li	•
Reach CROS SR	-
Reach CROS iX-CIC Li	-

Programming	
Noahlink Wireless	•
ConnexxAir	-
HI-PRO	-

• available - not available