# REACH SR




- Enabling users to be fully involved in group conversations with Multi-Voice Focus
- Hands-free phone calls via Headset Mode with Tap Control when connected via MFi or LE Audio
- Portable charger with built-in power bank
- Rechargeable solution for a full days use:
  - SR delivers up to 17h runtime (incl. 5h streaming)
- Lifeproof to withstand everyday challenges
- ATEX certified
- AI with Rexton Assist helps you to seamlessly customize your hearing aids

Colors

Black | Black Gloss


Cosmic Blue | Rose Gold

Snow White | Rose Gold

Snow White | Snow White Gloss


Black | Silver

### S-Receiver





- 46 dB / 110 dB SPL (2 ccm coupler)
- 56 dB / 120 dB SPL (Ear simulator)

### M-Receiver




- 60 dB / 119 dB SPL (2 ccm coupler)
- 70 dB / 129 dB SPL (Ear simulator)

### P-Receiver




- 65 dB / 122 dB SPL (2 ccm coupler)
- 75 dB / 131 dB SPL (Ear simulator)



| Features / Performance Level | 80 | 60 | 40 | 30 | 20 |
|---|---|---|---|---|---|
| Channels/Controls/Programs | 48/20/6 | 36/18/6 | 32/16/6 | 24/12/4 | 16/8/4 |
| LED | - | - | - | - | - |
| T-coil | - | - | - | - | - |
| Bluetooth | • | • | • | • | • |
| Tinnitus Sound function | • | • | • | • | - |
| Tinnitus Notch function | • | • | • | • | - |
| Wireless sync | • | • | • | • | • |
| IP68 | • | • | • | • | • |
| XPhone | • | • | • | • | • |
| **Technology Highlights** | | | | | |
| Multi-Voice Focus | • | • | • | - | - |
| *Focus Beams* | 4 | 2 | 1 | - | - |
| Voice Stabilizer | • | • | • | • | • |
| Speech Preservation Technology | • | • | • | • | • |
| My Voice 2.0 | • | • | • | - | - |
| Motion sensor for scene analyzis | • | • | • | - | - |
| **Directionality** | | | | | |
| Auto Dir | • | • | • | • | • |
| Directional Adaptive | • | • | • | • | - |
| iOmni Sound Locator | • | • | • | • | • |
| Auto iFocus 360 | • | - | - | - | - |
| iFocus 360 | • | • | - | - | - |
| Auto Focus 360 | - | • | - | - | - |
| Auto Stereo iLock | • | • | • | - | - |
| OneMic Directionality | - | - | - | - | - |
| **Noise Reduction** | | | | | |
| Noise Management | • | • | • | • | • |
| Sound Smoothing | • | • | • | • | • |
| Wind Noise Reduction | • | • | • | • | • |
| Auto Echo Reducer | • | - | - | - | - |
| Reverb Reducer | • | • | - | - | - |
| **Sound Quality & Control** | | | | | |
| Bandwidth extension | • | - | - | - | - |
| Bandwidth compression | • | • | • | • | • |
| Feedback preventer | • | • | • | • | • |
| Music Enhancer (presets) | 3 | 3 | 1 | 1 | - |
| Automatic Acclimatization | • | • | • | • | • |
| **Audio Streaming** | | | | | |
| Bluetooth Classic | - | - | - | - | - |
| LE Audio | • | • | • | • | • |
| MFi | • | • | • | • | • |
| ASHA | • | • | • | • | • |
| **Headset Mode** | | | | | |
| Bluetooth Classic | - | - | - | - | - |
| LE Audio | • | • | • | • | • |
| MFi | • | • | • | • | • |
| **TeleCare** | | | | | |
| Full live remote tuning | • | • | • | • | • |
| Basic portal remote tuning | • | • | • | • | • |

| Accessories & Apps | |
|---|---|
| Remote Control | • |
| TV Sound | • |
| Smart Mic | • |
| Rexton App | • |

| Chargers | |
|---|---|
| Charging Station | - |
| Travel Charger | • |
| Charging+ Station | • |
| Universal Cable Charger | - |

| CROS | |
|---|---|
| Reach CROS R-Li T BC | - |
| Reach CROS R-Li | - |
| Reach CROS SR | • |
| Reach CROS iX-CIC Li | - |

| Programming | |
|---|---|
| Noahlink Wireless | • |
| ConnexxAir | - |
| HI-PRO | - |

• available - not available